


GaBP RING AUTOLOCK™ SYSTEM

Banded Roux-en-Y Gastric Bypass Surgery




The silicone rubber elastomeric GaBP RING is designed for internal use in the field of Bariatric obesity surgery. These bands are of a soft consistency that provide for patient comfort and safety. The bands are radiopaque and made from implant grade silicone rubber. The ends of the bands have a plastic one-way locking mechanism that can be applied manually or laparoscopically. A flexible introducer is included (preattached), and is used to guide the GaBP Ring through the lesser omentum and around the gastric pouch.

The band assembly is reinforced with a monofilament core for strength and to help control food intake. The Gastrostomy Site Marker is reinforced with polyester to aid in suture retention.

*Optional Gastrostomy Tube not included in Kit.

GaBP RING AUTOLOCK™ SYSTEM



ORDERING INFORMATION

Cat. No.	Description	Size
FPO48005-55	GaBP RING? Autolock Kit	5.5cm
FPO48005-60	GaBP RING? Autolock Kit	6.0cm
FPO48005-65	GaBP RING? Autolock Kit	6.5cm
FPO48005-70	GaBP RING? Autolock Kit	7.0cm
FPO48014-00	Gastrostomy Site Marker?	6.0cm

The Bariathec GaBP RING Autolock System is undergoing approval for sale in the USA.



Caution: This device is restricted for sale by or on the order of a physician.

***USA and International Patents Pending.**

Distributed by:
 Bariathec Corp.
 PO Box 4257
 Palos Verdes Peninsula, CA 90274
 Ph: 562.809.1020, Fax 562.402.2374
www.Bariathec.com / n_fobi@bariatec.com



EUROPEAN REPRESENTATIVE: OBELIS S.A.
Corporate Address: Bd. Général Wahis, 53 1030
 Brussels, Belgium **Registered Address:** Avenue
 de Tervuren, 34, bte 44 1040 Brussels, Belgium